




### 1. Top Page

① EmerSim (Emergency care training system by Simulated vital sign monitor) is a training system with a vital monitor simulator. Below is a picture showing [EmerSim] application top page and menu settings.



\* Displaying vital values are listed below:  
HR(Heart Rate), SpO2(Peripheral Oxygen Saturation), BP(Blood Pressure), RESP(Respiratory Rate)

\* Caution: To use [EmerSim] application, Wi-Fi® setup is required.  
Refer to "EmerSim - Quick Start Guide" for Wi-Fi® setup procedure.

### 2. Network display



Displaying currently connected network name.

### 3. Sound on/off

Press speaker icon to turn sound on or off.  
Sound enabled example:




Sound disabled example:



### 4. Monitor

#### 4-1. Displaying vital values

By pressing "Monitor" button on application Top Page, vital values will be displayed on the monitor screen.



\* For operating vital values shown on the monitor, refer to "5. Controller".

#### 4-2. Returning to Top Page


To return to Top Page from Monitor screen, press on "Back".



### 5. Controller


#### 5-1. Scenario Select

1) From Scenario List on Application Top Page, select Scenario to display on the Monitor.



(\* above is an image of example scenario list)

2) Press "Controller" button.



3) Operational screen for the Monitor will be displayed.




■ **5-2. Switching Scene**

1) On top-right of the page, press "SCENES" button to display scenes list. On the list, scenario description and vital values for each scene will be displayed.



2) By selecting the scene, vital values will be updated to predefined values.



■ **5-3. Operating vital values**

Vital values are operated by a slider, numeric keypad and "+/-" buttons.

■ **5-3-1. Operating with slider**



By moving slider to the left, vital value will decrease,  
by moving slider to the right, vital value will increase.

■ **5-3-2. Operating with numeric keypad**



By pressing on a vital value number, numeric key will be displayed.  
After inputting value, press "OK" to update the vital value.

■ **5-3-3. Operating with "+ / -" buttons**




By pressing "+" button, vital value will increase.  
By pressing "-" button, vital value will decrease.


■ **5-4. Operating Waveform**

1) Refer to picture below and press around the area to change waveform.

(Example image on PC)




(Example image for smartphone)




2) Waveform list will be shown. Select waveform from the displayed list.

(Example image for PC)




(Example image for smartphone)



- Available waveforms for Heart Rate
- 1) Normal - NSR (Normal Sinus Rhythm)
  - 2) 1°block (First-degree atrioventricular block)
  - 3) 2°block (Mobitz I) (Second-degree atrioventricular block)
  - 4) 2°block (Mobitz II) (Second-degree atrioventricular block)
  - 5) 3°block (Third-degree atrioventricular block)
  - 6) Asystole
  - 7) Atrial fibrillation
  - 8) Atrial flutter
  - 9) Junctional
  - 10) LBBB (Left bundle branch block)
  - 11) RBBB (Right bundle branch block)
  - 12) ST elevation
  - 13) SVT (Supraventricular Tachycardia)
  - 14) Torsades de pointes
  - 15) VF (Ventricular fibrillation)
  - 16) VT (Ventricular tachycardia)

- Available waveforms for SpO2
- 1) Normal
  - 2) Poor perfusion


3) The waveform will be changed. (Example : from Normal-NSR to VF)



■ 5-5. Switching to shown or hidden.

In monitor's initial state all items are hidden.

1) Select an item to be shown on monitor. In picture below "HR" display status is being switched.



2) By switching to show state, value and waveform will be displayed on the monitor screen.

By switching to hide state, value and waveform will be hidden on the monitor screen.




On the monitor's screen:  
HR value and waveform are displayed



On the monitor's screen:  
HR value and waveform are hidden

■ 5-6. Returning to Top Page


To return to Top Page from Controller screen, press on "Back".




6. Scenario and Scene settings

■ 6-1. Registering Scenarios and Scenes.


1) On application's Top Page, press on "Settings" -> "Scenario Settings"




2) In Scenario column press "+" button.




3) Input "Scenario Name", "Description" and press "SAVE".  
As Scenario saves, "Insert Completed" will be displayed.



4) To add scene, select Scenario from scenario column and Press "+" "+" in Scene column.




5) Input "Scene Name", "Description", each vital value and press "SAVE".  
As Scene saves, " Insert Completed " will be displayed.






6-2. Editing/Removing Scenarios and Scenes.

1) Select an item from Scenario or Scene column and press "Edit" button.



2) Scenario or Scene editing screen will be displayed.




3) To edit, input "Scene Name", "Description", each Vital value and press "save". As Scenario or Scene updates, "Update Complete" will be displayed.

4) To remove, press on "Remove" button.

[ Troubleshooting ]


Can't press "Remove" button.  
When Scenario includes Scenes, it's unremovable.  
Remove all included Scenes in the Scenario to remove it.



7. Editing Common Settings

7-1. Changing color of vital items and waveform.


1) On application Top Page, press on "Settings" -> "Common Settings"




2) Press on "Font Color" for target Vital item. In the example below, RESP font color is changed from white to red.

Measurement Item	Notation	Font Color	Alert Settings
1. Blood Pressure	BP	yellow	Upper Limit: 180, Lower Limit: 40, <input checked="" type="checkbox"/> Blinking Alert, Alarm Sound: Alarm 1
2. Heart Rate	HR	green	Upper Limit: 150, Lower Limit: 30, <input checked="" type="checkbox"/> Blinking Alert, Alarm Sound: Alarm 1
3. Peripheral Oxygen Saturation	SpO2	blue	Upper Limit: 94, Lower Limit: 94, <input checked="" type="checkbox"/> Blinking Alert, Alarm Sound: Alarm 1
4. Respiratory Rate	RESP	white	Upper Limit: 40, Lower Limit: 11, <input checked="" type="checkbox"/> Blinking Alert, Alarm Sound: Alarm 1

3) Select color from the list (as "red" in the example) and press "SAVE" button. "Update Complete" will be displayed.




4) Vital value font color and/or waveform color will be updated on Monitor and Controller screens.



7-2. Setting Vitals Alarm (blink / warning sound)

1) In Common Settings screen, edit "Alert Settings" values for specific vital. Update Blink upper/lower limit and check checkbox to enable it.

Measurement Item	Notation	Font Color	Alert Settings
1. 血圧	BP	yellow	Upper Limit: 180, Lower Limit: 40, <input checked="" type="checkbox"/> Blinking Alert, Alarm Sound: Alarm 1
2. 心拍数	HR	green	Upper Limit: 150, Lower Limit: 30, <input checked="" type="checkbox"/> Blinking Alert, Alarm Sound: Alarm 1
3. 動脈血酸素飽和度	SpO2	blue	Upper Limit: 94, Lower Limit: 94, <input checked="" type="checkbox"/> Blinking Alert, Alarm Sound: Alarm 1
4. 呼吸数	RESP	white	Upper Limit: 40, Lower Limit: 11, <input checked="" type="checkbox"/> Blinking Alert, Alarm Sound: Alarm 1




Without settings upper limit, even by checking blink, it won't work upon reaching the limit. Same with lower limit. Without checking Blink checkbox, values won't blink upon reaching set upper/lower limit.

2) Press "SAVE" button to complete update.





3) Upon reaching lower or upper limit in Monitor or Controller screen, target Vital value will blink.



**7-3. Showing/hiding logo on Monitor screen**


**Logo is displayed on Monitor screen by default**

1) In "Common Settings", find "Monitor Screen Logo" [Other Settings] and press on select arrow. In the example below, Logo is being switched off.




**7-4. Switching screen output (HDMI) settings**

1) On application's Top Page, press on "Settings" -> "Screen Output Settings"




2) It's possible to change HDMI display settings for [EmerSim] device depending on TV/PC monitor. Select connecting screen type ("Connect to TV") or ("Connect to PC monitor"), and select according "Resolution, Refresh Rate".

\* Check your display manual for "resolution" and "refresh rate".




3) In the example below, "PC monitor connect" is selected with "Resolution, Refresh Rate" of [20. 1024x768 120Hz].



4) After selecting "Resolution, Refresh Rate" press "APPLY SETTINGS" button. A confirm message saying ["To apply settings, the device will be restarted. During restart Wi-Fi connection will be lost. Please reconnect to the device after restart is complete."] will be displayed. Press "OK" button.



5) [EmerSim] device reboots and Wi-Fi connection will be lost as on example below. Refer to "EmerSim - Quick Start Guide" [3-2. Preparing your PC, smartphone] for Wi-Fi® setup procedure.



**Network Error**  
Disconnected from [qtrain8001]. Check Wi-Fi connection to the device and ensure the cable is properly plugged in.

6) As Wi-Fi® connects, Top Page will be displayed with set resolution and refresh rate".






## 8. Optional Functions

### 8-1. Simulation Sound On/Off/Mute

- 1) Along with monitor sound, there are playable simulation sounds. By pressing on "Bell Icon" on the top, playlist will be displayed.



- 2) In the example below, "phone" simulation sound is being played.



- 3) To stop simulation sound, press ■ "Stop" button.
- 4) To mute simulation sound, press 🚫 "Mute" button.

### 8-2. Stopwatch function

- 1) Press "Start" to begin displaying elapsed time.
- 2) Press "Stop" to stop displaying elapsed time.
- 3) Press "Reset" to reset elapsed time.



## 9. Troubleshooting

Please check content below for troubleshooting.

### 9-1. [EmerSim] device is not working.

Is the device turned on?

Refer to "EmerSim - Quick Start Guide" [3-1. Preparing [EmerSim] device] and check if device LED is red.

Is the device overheated?

[EmerSim] device might not function properly under high temperatures. Unplug the device and let it cool down before turning it on again. Do not place the device near heat generating sources or in direct sun.

Is Wi-Fi® connected?

Refer to "[9-2. Wi-Fi® Connection Issues](#)".

### 9-2. Wi-Fi® Connection Issues

Is Wi-Fi® signal strong enough?

Being away from the [EmerSim] device or having wall in between might cause connection issues. Try connecting while being closer to the device.

Is Wi-Fi® setup complete?

Refer to "EmerSim - Quick Start Guide" [3-2. Preparing your PC, smartphone] and verify settings.

### 9-3. Vital items / waveforms are not displayed on the monitor

In Controller set each vital item display setting to "ON". Each vital item has a monitor icon. Please click on the icon to switch between show and hide mode.

By default, vital items are hidden on the monitor.

Please set each vitals visibility according to the training procedure.

Refer to "[5-5. Switching between show and hide state.](#)"

### 9-4. Sound is not working

Does connected display support speakers?  
Check if speakers function properly.

Set sound to On.

Refer to "[3. Sound on/off](#)"

### 9-5. HDMI output is not working

Is device resolution and refresh rate set according to connected TV or PC monitor?

Refer to "[7-4. Switching screen output \(HDMI\) settings](#)"

## 10. Usage Precautions and Warranty Period

Please read articles below in "EmerSim - Quick Start Guide" and accept usage agreement.

- [6. Usage Precautions]
- [7. Radio Wave Precautions]
- [8. Power Source Precautions]
- [9. Main Specifications]
- [10. Warranty Period]

Upon usage, user accepts agreement to above mentioned.

"Seller warrants to Buyer that the Product will materially conform to the specifications set forth in Exhibit.

EXCEPT FOR THE WARRANTY SET FORTH ABOVE, SELLER MAKES NO WARRANTY WHATSOEVER WITH RESPECT TO THE PRODUCT, INCLUDING ANY (A) WARRANTY OF MERCHANTABILITY; OR (B) WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE; WHETHER ARISING BY LAW, COURSE OF DEALING, COURSE OF PERFORMANCE, USAGE OF TRADE OR OTHERWISE."

## Inquiries

Developer/Distributor  PENGUIN SYSTEM® Phone Number +81-29-893-2275  
 ペンギン システム 株式会社 URL <https://www.penguins.co.jp/en>  
 E-mail [qtr@penguins.co.jp](mailto:qtr@penguins.co.jp)



publication number: Ver.1.0.5-r4

